


PHT (15-25) 3.5L 60 Hour - Packing Instructions



Assembly Instructions

Prepare two BIOCHILL 1515s, and four BIOCHILL 2020s by conditioning them at +20°C for a minimum of 24 hours.

Once the BIOCHILLS are fully conditioned, place one ambient BIOCHILL 1515 into the centre of the EPS 10 base, and add the filled payload carton directly on top of the BIOCHILL.

Next add an ambient BIOCHILL 2020 around each side of the payload carton. Next stack the remaining ambient BIOCHILL 1515 on top of the payload carton, as pictured.

Finally apply the EPS 10 box lid, then close and seal the outer carton. The system is now ready to be shipped.